FIN SYSTEM

COTILI ACETABOLARI
ACETABULAR CUPS
COPAS ACETABULARES

FIN SYSTEM

COTILI ACETABOLARI
ACETABULAR CUPS
COPAS ACETABULARES

FIN SYSTEM

IT

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FIN SYSTEM

EN

The FIN SYSTEM acetabular system represents a highly versatile solution with a wide choice of liners and heads, also of large diameters. The cups are also united by the same surgical procedure and instrument kit.

The system includes 4 cups with specific characteristics and indications:

- FIN II: uncemented primary cup, characterized by three cranial fins; made of Titanium alloy Ti6Al4V grade 5 ELI (ISO 5832-3).
- FIN MB: uncemented primary cup; made of Titanium alloy Ti6Al4V grade 5 ELI (ISO 5832-3).
- FIN DMD: uncemented primary cup, characterized by a lattice structure that replicates that of a diamond; made of Titanium Ti6Al4V ELI (ASTM F3001; ISO 5832-3) by additive manufacturing.
- FIN DMD MULTIHOLE: uncemented revision cup, characterized by a lattice structure that replicates that of a diamond; equipped with a wide range of wedges; made of Titanium Ti6Al4V ELI (ASTM F3001; ISO 5832-3) by additive manufacturing.

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Liners

EN

The compatible liners come in four versions:

- Ultra-heavy molecular weight polyethylene (UHMWPE - ISO 5834-2) liners with or without antiluxation shoulder; int. diameter 28mm.
- XLPE cross-linked polyethylene (ISO 5834-2) liners with or without antiluxation shoulder; int. diameter 28mm, 32mm and 36mm.
- XLPE liners with vitamin E (ISO 5834-2) with or without antiluxation shoulder; int. diameter 28mm, 32mm and 36mm.
- BioloX® DELTA ceramic liners (ISO 6474-1,-2); int. diameter 28mm, 32mm, 36mm, 40mm*.

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
FIN II

EN

The FIN II cup features three peripheral fins to improve primary stability and to prevent rotation.

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
Three holes in the proximal zone of the cup for optional application of 6.5 mm screws. The specially designed screw housings provide axial deviation of up to 15°

FIN II: schrew hole plugs provided for all sizes except sz. 44
FIN MB: provided without schrew hole plugs

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A polar hole for direct visual inspection of the cup position with respect to the bottom of the acetabulum

FIN II: Apex pug provided for all sizes except sz. 44
FIN MB: provided without apex pug

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Cranially Placed at 70° from each other

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FIN MB

EN

The FIN MB features the same morphological characteristics of the FIN II cup except for the three cranial fins. Though it maintains primary and secondary stability thanks to its geometry and the porous coating.

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Geometry

Hemispherical, allowing global contact between the component and the bone acetabulum

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
FIN DMD

EN


Uncemented primary cup, it is characterized by a lattice structure that replicates that of a diamond.

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Tecnologia Additive Manufacturing

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Additive Manufacturing Technology

Additive Manufacturing technology makes it possible to produce a cup in a single process. It is characterized by a lattice structure that replicates that of a **DIAMOND**.

The size and morphology of the cells provide an ideal habitat for the osteoconduction process, thus accelerating osseointegration.

Tecnología Additive Manufacturing

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Holes

The holes are intended for the use of FIN system cancellous screws; diameter 6.5 mm.

FIN DMD: 3 holes

FIN DMD MULTIHOLE: 10 holes in sizes 46, 48 and 50 and 12 holes in sizes from 52 to 72

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
FIN DMD MULTIHOLE

EN

Uncemented revision cup.
The cup is equipped with a wide range of wedges allowing the surgeon to deal with most acetabular revision cases, even in the presence of cavitory and segmental bone defects.

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Material


FIN DMD and **FIN DMD MULTIHOLE**: Titanium Ti6Al4V ELI (ASTM F3001; ISO 5832-3) by Additive Manufacturing technology

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WEDGES

EN

The acetabular wedge system for the FIN DMD MULTIHOLE cup was designed for surgical operations of the hip that require extra support for the acetabular element in cases where the bone stock is insufficient.

- Intraoperative flexibility
- Designed for non-cemented fixation on the bone/wedge interface and for cemented fixation on the cup/wedge interface
- Three thicknesses (15, 20 and 25mm)

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
UHMWPE LINER

REF	SIZE	DIAM
110364242	42mm	28mm
110364244	44/46mm	28mm
110364248	48/50mm	28mm
110364252	52/54mm	28mm
110364256	56/58mm	28mm
110364260	60÷72mm	28mm

XLPE ANTILUXATION SHOULDER LINER

REF	SIZE	DIAM
110361142	42mm	28mm
110361144	44/46mm	28mm
110361248	48/50mm	32mm
110361252*	52/54mm*	32mm*
110361352	52/54mm	36mm
110361356	56/58mm	36mm
110361360	60÷72mm	36mm

XLPE + VIT E ANTILUXATION SHOULDER LINER

REF	SIZE	DIAM
110363342	42mm	28mm
110363344	44/46mm	28mm
110363348	48/50mm	32mm
110363352	52/54mm	36mm
110363356	56/58mm	36mm
110363360	60÷72mm	36mm

UHMWPE ANTILUXATION SHOULDER LINER

REF	SIZE	DIAM
110364342	42mm	28mm
110364344	44/46mm	28mm
110364348	48/50mm	28mm
110364352	52/54mm	28mm
110364356	56/58mm	28mm
110364360	60÷72mm	28mm

CERAMIC BIOLOX® DELTA LINER

REF	SIZE	DIAM
110367044	44/46mm	28mm
110367048	48/50mm	32mm
110367152	52/54mm	36mm
110367156	56/58mm	36mm
110367176*	56/58mm*	40mm*
110367160	60÷72mm	36mm
110367180*	60÷72mm*	40mm*

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