# “M.R.P.”

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**COMPONENTE  
DISTALE**  
*DISTAL  
COMPONENT*

**COMPONENTE  
PROSSIMALE**  
*PROXIMAL  
COMPONENT*

## **A MODULAR SYSTEM TO ASSEMBLE INTRAOPERATIVELY**

Limb salvage surgery is today a reliable option in the management of bone tumors. During the past 25 years a large experience in proximal femur resection and prosthetic replacement has been developed by many centers throughout the world. (1-5).

Goals of proximal femur replacement are preservation of limb length, restoration of appropriate hip center, and achievement of satisfactory stability throughout a functional range of motion. The Modular Resection Prosthesis (M.R.P.) system combines simplicity and versatility providing a wide range of reconstructive options. A 135° neck-shaft angle improves the abduction lever arm optimizing hip stability. Modularity allows intraoperative assembly customizing the implant to each specific need.

## **INDICATIONS**

Current indications to proximal femur replacement for neoplastic conditions include:

- primary malignant bone tumors not involving major neurovascular structures;
- metastatic cancer;
- benign bone tumors presenting with extensive soft tissue expansion (stage 3) or pathologic fracture.

## **CHARACTERISTICS**

The system consists of three modular basic elements:

- the proximal component (135° neck - shaft angle);
- the distal component;
- the locking screw;

## **MATERIALS**

- Titanium Ti - 6Al - 4V - grade 5 - ELI ISO 5832/3

## **MODULARITY**

The modularity of each basic element of the system is as follows:

- Proximal component (75 - 100 - 125 - 150 - 175 - 200 mm) length;
- Distal component (140 - 160 - 180 mm) length and (12 - 13 - 14 mm) diameter;
- Locking screw (75 - 100 - 125 - 150 - 175 - 200 mm)

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Male, age 71, grade 2 chondrosarcoma of the proximal femur.

## TECNICA OPERATORIA - OPERATING TECHNIQUE

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**A)** The patient is usually placed in a lateral decubitus position and the resection then carried out as dictated by each specific need. The resection specimen length is measured and compared to preoperative planning. Limb length restoration or 1-1,5 cm overlengthening should be accomplished according to the amount and quality of remaining soft tissue, in order to obtain satisfactory implant stability.

**B)** The shortest body allowing length restoration should be selected in order to preserve bone. Stem length and diameter are selected according to preoperative planning and intraoperative findings. In most circumstances we prefer bipolar to total joint replacement. Proximal and distal component are assembled and then locked by the locking screw using the blocking wrench to accomplish secure fixation. Once the appropriate bipolar has been chosen, a trial can be carried out using a standard head (0 mm), checking leg length and stability. Additional 0,5 - 1 cm of bone can be sacrificed to obtain the desired length and soft tissue tension.

**C)** The medullary canal is then prepared and the implant cemented according to the surgeon's preferred technique. A new trial should be carried out to check again length and stability; the most satisfactory modular head is then selected and the real impacted.

**D)** The hip is reduced and closure undertaken according to the surgeon's preference. An effort should always be made to maximize soft tissue tension by satisfactorily reattaching the remaining abductors to the fascia and/or prosthesis.

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GRUPPO BIOIMPIANTI S.r.l.