SMR

STELO MODULARE DA REVISIONE
MODULAR HIP REVISION STEM
VASTAGO MODULAR DE REVISION

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SMR MODULAR HIP REVISION STEM

EN

The SMR hip revision stem is a modular system designed for prosthesis replacement in cases of serious bone loss ("Società italiana riptotesizzazione-GIR" grades II, III and IV and Paprosky grades II and III). With over 20 years of clinical history, the femoral implant consists of three elements: a proximal component, and distal component and locking screw.

The SMR stem allows customized reconstructions thanks to the wide range of modular combinations. The distal component of the SMR revision stem is also used in coupling with the proximal component of the MRP modular resection prosthesis. The surgeon has thus the possibility of using cemented or uncemented distal components for both SMR and MRP prosthesis, with different indications.

Indications

EN

In this surgical technique are explained the guidelines for SMR Modular Revision stem implant. The SMR stem is indicated in: revision surgery of femoral components, following failure of primary cemented or un-cemented prosthesis with proximal bone loss, enlargement of the medullary canal and thinning of the cortical bone of the proximal femur; revision of a mobilized femoral component, after sub prosthetic or periprosthetic fracture; femoral reconstruction in cases of alterations of bone morphology after fractures or osteotomies; this stem is indicated for cementless implants.

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Contraindications

EN

The hip joint surgery is absolutely contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis. It is relatively contraindicated in case of: osteoporosis; patient uncooperative or suffering from neurological disorders, unable to follow directions; systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; osteomalacia; active infection or suspected latent infection in the hip joint; distant focus of infection that could spread to the implant site; persistent acute or chronic osteomyelitis; vascular insufficiency, muscular atrophy, neuromuscular diseases; incomplete or insufficient presence of soft tissue around the knee joint; obesity; inadequate bone stock for the prostheses support or fixation; skeletal immaturity local or disseminated neoplastic diseases; incorrigible severe deformities.

Materials

EN

Titanium alloy Ti 6Al 4V Grade 5 ELI, which ensures excellent bio-compatibility and mechanical strength. The surface finish is obtained by abrasion with corundum to increase the contact surface and thereby improve prosthesis-bone anchorage.

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PROXIMAL COMPONENT

EN

- The modular proximal component solves the problem of restoring the physiological length of the limb and stabilizing the joint thanks to its four lengths 40, 50, 60 and 70 mm and two CCD angles 125° and 135°.
- The mirror polished neck area has been studied to decrease polyethylene wear during flexion. In addition, the mirror finish eliminates the risk of fatigue cracks.

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
DISTAL SCREWS

The 160, 200 and 240 mm distal components have three holes, positioned radially in the distal portion at 120° intervals, to take 4.5 mm cortical screws in the rare cases they are necessary.

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FILLING PROXIMAL COMPONENT

EN

Given the implant's modularity and versatility, the SMR system offers several advantages. It also makes it possible to opt for mounting filling proximal components on the distal components, which is useful in all cases of:

- Removal of massive quantities of bone following removal of the previously placed implant;
- Difficulty in removing the stem making large resections necessary;
- Formation of bone defects due to regressed sepsis or osteolysis caused by debridement of prosthetic materials.

The filling proximal component is also characterized by the presence of a further lateral shoulder with three holes for passing wires or cables (up to 3 mm diameter) for cerclage when resolving cases of fractures/detachments of the greater or lesser trochanter.

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DISTAL COMPONENT

EN

- The distal component is conical and has longitudinal fins, this prosthesis design guarantees axial stability and resistance to torsional forces at the same time.
- The modular design has four lengths (100, 160, 200 and 240 mm) and ten diameters (14, 15, 16, 17, 18, 19, 20, 22 and 24 mm).
- The 200 e 240 mm distal components have a 3° anatomical pro-curvedness to avoid point effects on the cortical.
- The distal portion has three holes placed radially at 120° for cortical screws to be used in the rare cases they are needed.

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Pre-operative planning

EN

Thorough pre-operative planning enables the surgeon to adopt a correct approach to the operation and foresee certain situations that might arise.

The objective is to choose the diameter and length of the distal component and size of the proximal component. Pre-operative choice of components is indicative, as the definitive size will be determined intra-operatively.

Pre-operative planning requires an anterior-posterior X-ray of the pelvis and pre-operative templates of the SMR stem (provided in the instrument kit). The templates are 15% enlargements.

Align the template with the anatomical axis of the femur to determine the size of the distal component, making sure that it's well anchored to good quality cortical bone.

Any anatomical differences may be corrected using the femoral heads and proximal components, which come in different lengths.

Also determine the exact position of the longitudinal osteotomy, where possible.

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Reaming of distal component

EN

Having completed the removal of the previous implant, proceed to reaming with the conical reamers (Ref. 120540514 ÷ 120540524*, 120540714 ÷ 120540724*). (FIG. 2a, FIG. 2b)

Reaming may be done manually with the "T" handle (Ref. 110381521) or with a motor, using the Jacobs Adaptor (Ref. 110381525).

The instrument kit has reamers with diameters from 14 mm to 20 mm, in 1 mm steps, and 22 mm and 24 mm.

Start with a reamer that's one or two millimetres smaller than the measurement taken in pre-operative planning.

The reamers must be inserted into the femoral canal till the depth indicator is aligned with the apex of the greater trochanter.

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EN

The shaft of the reamers L.160 (Ref. 120540514 ÷ 120540524*) is graduated and divided into three groups of numbers, from 40 to 70, corresponding to the four available sizes of the proximal components, while the three groups refer to the three lengths of the distal stem (160, 200, 240). (FIG. 2b)

For example, if the reamer stops at 60 in the first groups of numbers (the one closest to the tip of the reamer), a 160 mm long stem with a 60 mm proximal component will be used. The shaft of the reamers L.100 (Ref. 120540714 ÷ 120540724) is graduated and divided from 40 to 70, corresponding to the four available sizes of the proximal components. (FIG. 2a) Work with conical reamers of successive sizes till achieving optimum stability and resistance and coming into contact with the cortical.

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Insertion of distal component

EN

Use of the reamers yields information on the diameter and correct depth of the distal component in the cortical.

Mount the distal component inserter (Ref. 120540410) on the chosen component and proceed to feed it into the previously prepared canal.

The long distal components (200 mm and 240 mm) must aligned with the femoral pro-curvedure (3°).

Proceed to insertion with the distal component inserter advancing gradually and checking the depth, as far as the maximum, to obtain good torsional stability and avoiding any unexpected penetration under load.

(FIG.3)

Indicatively, the final position of the component will correspond to the mark identified on the reamer in the previous phase.

Any difference may be corrected using the proximal components and femoral heads, which come in different lengths.

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
FIG. 4a

Preparation of the bone for the proximal component

EN

Remove the distal component inserter (Ref. 120540410) from the distal component, screw the guide shaft for tubular reamer (Ref. 120540415) onto the distal component. (FIG.4a) When carrying out this operation, make sure it's done in the right direction, as a mistake could irreparably damage the distal component and compromise all the successive phases. Now ream the proximal portion with the appropriate tubular reamer (Ref. 120540400) to prepare a housing for the proximal component. (FIG.4b)

Attention: proceed with prudence in the bone reaming phase to avoid damaging any soft tissue in the field of action.

Attention: the reamer must rest on the apex of the Guide in order to avoid difficulties with insertion in the subsequent housing of the proximal components. To be certain of carrying out the procedure correctly, make sure there is no residue inside the reamer that might hinder its correct penetration in the bone. There is an inspection hole on the reamers for checking there are no residues.

FIG. 4b

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Insertion of trial proximal component

EN

Insert the trial proximal component (Ref. 12540140 ÷ 1205401701; 12540240 ÷ 120540270) and lock it with the appropriate screw (Ref. 120540040 ÷ 120540070). Do not tighten excessively. (FIG.5a)

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Carry out the trial reduction and assess the stability of the hip, the length of the limb and the breadth of movement. Also check for the desired anteversion, marking the bone by the reference line marked laterally on the trial proximal component.

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Insertion of the definitive proximal component

EN

Remove the trial head and trial proximal component using the revision proximal component extractor (Ref. 120540420).

Proceed with prudence in the extraction phase and make sure the correct direction is chosen, as a mistake could compromise all the successive phases.

Take the definitive proximal component (having the same size as the trial one) and use the revision proximal component impactor (Ref. 120540411) to insert it in the Morse taper of the previously cleaned distal component. (FIG.6)

Be careful to observe the correct anteversion by aligning the reference line on the proximal component with the marking previously made on the bone. Definitively impact the proximal component. The filling proximal component does not require any steps beyond those indicated in the proximal component technique.

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Blocking the proximal component

EN

Insert the locking screw (of the appropriate size for the proximal component used) and tighten it with the screwdriver for revision/resection prostheses (Ref. 120540405). When tightening the screw, use the non-rotation device (Ref. 120540412), which fits into the Morse taper, to prevent rotation during assembly. (FIG.7)

The locking screw ensures additional anchoring between the distal and proximal components, as well as holding the Morse taper.

Insert the definitive head (corresponding to the trial used) by impacting it with the appropriate adaptor for head impactor 36mm (Ref. 110380812) mounted on the head impactor (Ref. 110380800) and reduce the prosthesis.

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Extraction of SMR proximal component

EN

Should it be necessary to extract the proximal component of the SMR stem, remove the locking screw, screw the extractor for revision proximal component (Ref. 120540420) into the locking screw hole and turn until it's fully extracted. (FIG.8)

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Extraction of SMR distal component

EN

Should it be necessary to extract the distal component of the SMR Stem, use the instruments provided for the purpose. Screw the extension for distal component extraction (Ref. 120540421) Onto the distal component and use the slap hammer (Ref. 140030665) to extract it.

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