# JANUS CUP

CUPOLA BIARTICOLARE  
BIPOLAR CUP  
COPA BIARTICULARE

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## JANUS CUP

EN

The bipolar bicentric JANUS CUP ensures biomechanical advantages that monocentric bipolar cups do not have. When head and cotyle have different centres of rotation they spontaneously tend to align themselves under load along the longitudinal axis of the neck, which coincides with the line of action of the force weight.

The result is that in a bicentric bipolar the mobile cup will never slip out of line and will always stay well centred with respect to the line of load. This means less risk of dislocation, reduced distribution of load on the natural cotyle, less likelihood of seizure and, last but not least, enhanced radiography acceptability.

## Indications

EN

Bipolar JANUS Cup is indicated in hemiarthroplasty on base cervical fractures of the femoral neck and they articulate directly with the natural acetabulum.

The general guidelines are: non-inflammatory degenerative joint disease such as primary or secondary osteoarthritis; aseptic necrosis of the femoral head; rheumatoid arthritis; post-traumatic Arthritis; correction of functional deformity; outcomes of fractures of the femoral neck; outcomes of traumatic dislocations of the hip; failures of osteotomy; outcomes of arthrodesis.

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## Contraindications

EN

The hip joint surgery is specifically contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis. It is relatively contraindicated in case of: osteoporosis; patient being uncooperative or suffering from neurological disorders, unable to follow directions; systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; osteomalacia; active infection or suspected latent infection in the hip joint; distant focus of infection that could spread to the implant site; persistent acute or chronic osteomyelitis; vascular insufficiency, muscular atrophy, neuromuscular diseases; incomplete or insufficient presence of soft tissue around the hip joint; obesity; inadequate bone stock for the prostheses support or fixation; skeletal immaturity; local or disseminated neoplastic diseases; non treatable severe deformities.

## Materials

EN

**CUP:**  
Stainless steel to ISO 5832-1 standard.  
**INSERT:**  
Ultra high molecular weight polyethylene (UHMWPE) GUR 1020 to ISO 5834-2 standard.

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### INSERT

Ultra high molecular weight polyethylene (UHMWPE)  
GUR 1020 to ISO 5834-2 standard.

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
## JANUS CUP

EN

The tapered profile of the metal edge of the cup stops it coming into contact with the cotyle's articular cartilage, even in extreme positions. The "grating" effect is thus avoided and the cotyle wear typical of endoprosthesis is minimized.

## Heads

The Inserts of Janus Cup are compatible with heads of 28mm diameter.

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## Features

EN

- Less likelihood of seizure;
- Reduced load distribution on the natural cotyle;
- Enhanced radiography acceptability.

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$$M = F_s \times D$$

## Bicentric

In a BICENTRIC SYSTEM the centre of the head and cup align themselves along the axis of the forces being exerted.

Elimination of moment due to alignment of forces is obtained when  $D = 0$ . This guarantees the system greater stability.

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

D= Torque arm or distance between the two centers  
Fp= Force weight acting on the centre of the cup, downwards  
Fs= Supporting force exerted by the limb acting on the centre of the head, upwards  
m= Resulting torque

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### Pre-operative assessment

EN

A bilateral X-ray of the pelvis will help gauge the length of the limb, which has to be determined pre-operatively and then corrected.

The pre-operative templates of JANUS cups are provided in a 15% enlarged format.

Use of an X-ray reference may help to determine the enlargement of the patient's X-ray.

It is possible to obtain the measure of the head of the explanted femur, through the use of a caliper.

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### Trial implant assembly

EN

Once the size of the acetabulum has been determined, insert the trial head into the trial cup (Cod. 110450240 ÷ 110450258).

(FIG. 1)

Mount the assembly on the taper of the previously implanted stem and proceed with the trial reduction, to verify the correct choice of the cup size and the correct length of the head. (FIG. 2)

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#### Trial implant disassembly

EN

To disassemble the trial head from the trial cup, use the Trial head extractor (Ref. 110460101): insert the trial head extractor in the head and remove the head, holding the trial cup with one hand. (FIG. 3)

Continue the movement until the two trial components are detached.

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#### Definitive implant assembly

EN

To proceed with the assembly of the final implant:

1. Insert the definitive head into the polyethylene insert of the size corresponding to that of the definitive cup. (FIG. 4A)
2. Insert the insert-head assembly into the final cup and proceed to the implant. (FIG. 4B)

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Should it be necessary to disassemble the final implant, use the Insert extractor (Ref. 110460100), inserting the two pins into the holes in the polyethylene insert. (FIG. 5)

Tighten the extractor-pincer, closing it until the insert is removed from the metal cup.

If necessary, to extract the definitive head from the polyethylene insert, use the Trial head extractor (Ref. 110460101).

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


### Sterilization

EN

Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1,-2,-3 while the elements manufactured in UHMWPE are sterilized with ethylene oxide (ETO) according to the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioner.

### Identification and traceability

EN

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method. To guarantee external traceability even after implantation they are provided six labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

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## Warnings

EN

Never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant. In particular, never use stems and femoral heads produced by different manufacturers, because coupling tapers can be non-compatible. Gruppo Bioimpianti S.r.l. will not respond if any problems arise by a similar choice. The surgeon must be aware that a damaged or scratched neck can influence stem resistance and lead to a premature neck fracture.

1. Handle the product with due care to avoid accidental scratches, even during surgery. An inadequate handling of the device, pre-operative and / or intra-operative, or any damage can lead to corrosion, fatigue failure or excessive wear. Use clean gloves when handling the device.
2. Use screws or accessories made of suitable material: in the case of a system composed of several elements, always use the same kit elements: in case of acetabular cups made of titanium alloy (ISO 5832/3) use, if necessary, only screws of the same material.
3. Follow the surgical technique provided by Gruppo Bioimpianti.
4. A complete cleaning and removal of metal particles, bone cement, as well as any other debris resulting from surgery, it is necessary to minimize the wear of the articular surface.
5. Do not use the device if there is any damage.
6. Do not alter the system in any way.
7. Inadequate choice, positioning, alignment and fixation of one or more components could result in an incorrect distribution of the stresses, with a consequent reduction of life of the prosthetic implant.
8. An incorrect alignment or a plant not accurate can lead to a phenomenon of excessive wear and / or failure of the system itself.
9. All implantable devices for hip replacement are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use, which affects the fatigue strength and thus increases the risk of breakage.
10. Not implant a component in a patient if it had already been implanted, albeit temporarily, in another patient.
11. As regards the acetabular wedges of FIN DMD system fix the wedge/cup interface using the cementless fixation and fix the wedge/remaining acetabular bone using the cementless fixation by press-fit and 4.5mm screws, which must be passed through the holes present on the same wedge. Never use these screws between the two components (through cup and wedge): the wedge need to be fixed to the cup only using bone cement, without the use of the screws.
12. It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability.
13. Patient smoking may result in a delayed healing, non-healing and/or compromised stability.
14. It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components.

Implantable devices for hip replacement surgery provides the surgeons the necessary tools to reduce pain and restore the hip joint function. While this product is generally successful in attaining these goals, it cannot be expected to withstand the activity levels and function of a normal healthy joint.

## Precautions

EN

The instruments set, realized by Gruppo Bioimpianti, was designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades...) could determine negative events. Any fracture or breaking of instruments shall be reported. Surgical instruments are inevitably subject to wear with normal usage and as they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti. It is recommended that all instruments be regularly inspected to identify any possible damage. Do not leave in the surgical site any trial, packaging, instrument or its part.

## Possible adverse effects

EN

The possible adverse effects that may result from the use of a prosthetic system, are the classics that may result from a hip arthroplasty and include the following:

1. Osteolysis (progressive bone resorption). Osteolysis can be asymptomatic and therefore routine periodic radiographic examination is essential to prevent any serious future complication
2. Particles leading to increased wear rates necessitating early revision.
3. Early or late deep infection which may necessitate the removal of the prosthesis.
4. Allergic reactions.
5. Intra-operative bone fractures particularly in the presence of low bone density due to osteoporosis, bone defects resulting from previous surgeries, bone resorption, or insertion of the device.
6. Damage to blood vessels or hematoma
7. Dislocation, migration and / or subluxation of the prosthetic components due to trauma, incorrect attachment, loss of fixation, misalignment, improper positioning, bone resorption, soft tissue laxity, sudden movement and / or natural and / or excessive activity.
8. Peri-articular calcification or ossification, with or without impediment to joint mobility.
9. Undesirable variation of the length of the limb.
10. Reduced Range of Movement due to an incorrect choice of the components or a inadequate placement of the implant, or femoral impingement.
11. Fatigue fracture of the prosthetic components can occur as results of loss of fixation or incomplete fixation, strenuous activity, trauma, misalignment or excessive weight.
12. Corrosion at the interface between components.
13. Wear and / or deformation of the articular surfaces.
14. Pain.
15. Post-operative bone fracture post-operative and / or post-operative pain.
16. Audible sounds during motion.
17. Nerve damage.
18. Trochanteric avulsion or non-union as result of excess muscular tension, early weight bearing, or inadequate reattachment.
19. Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialization or muscle deficiencies.

## Information

EN

For any information or enquires about this publication or anything else, contact GRUPPO BIOIMPIANTI.

## Disclaimer

EN

This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons. This document does not constitute medical advice, it does not dispense medical recommendations and it does not convey any diagnostic or therapeutic information.

Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists, however Gruppo Bioimpianti excludes any liability for improper use of informations.

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