# FIN SYSTEM

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OPERATING TECHNIQUE

TÉCNICA QUIRÚRGICA

IT | EN | ES

*Enjoy Mobility*

# FIN SYSTEM

## FIN SYSTEM

IT

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## FIN SYSTEM

EN

The FIN II cup was designed in 1992 to offer a highly versatile solution with a wide choice of inserts and heads. The range has been enriched over the years to meet the specific needs of the surgeon. The system has grown and now includes 4 cups with specific features and indications. These cups are compatible with the system's inserts and share the same surgical procedure and instrument kit.

## Indications

EN

The FIN II, FIN MB and FIN DMD cups are indicated for cases of primary total arthroplasty to reduce or eliminate pain and/or improve hip functionality. The specific clinical indications are as follows: non-inflammatory degenerative articular pathologies, such as primary or secondary osteoarthritis; aseptic necrosis of the femur; rheumatoid arthritis; post-traumatic arthritis; functional deformity corrections; hip dysplasia; results of femur neck fracture; results of traumatic hip dislocation; failed osteotomy; results of arthrodesis. The FIN II, FIN MB and FIN DMD cups are recommended for uncemented use because of their porous 3-dimensional surface structure.

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The FIN DMD Multihole cup is indicated for total hip replacement revisions necessitated by the following conditions: mechanical failure of the primary prosthesis; instability of the primary prosthesis; periprosthetic fracture of hip; wear/osteolysis of the acetabular insert; aseptic mobilization; deep infection. The FIN DMD Multihole cup is recommended for uncemented use because of its porous 3-dimensional surface structure.

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## Contraindications

EN

The hip joint surgery is specifically contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis. It is relatively contraindicated in case of: osteoporosis; patient being uncooperative or suffering from neurological disorders, unable to follow directions; systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; osteomalacia; active infection or suspected latent infection in the hip joint; distant focus of infection that could spread to the implant site; persistent acute or chronic osteomyelitis; vascular insufficiency, muscular atrophy, neuromuscular diseases; incomplete or insufficient presence of soft tissue around the hip joint; obesity; inadequate bone stock for the prostheses support or fixation; skeletal immaturity; local or disseminated neoplastic diseases; non treatable severe deformities.

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## FIN II CUP

EN

Uncemented cup made of Ti6Al4V grade 5 ELI (ISO 5832-3), with a with titanium Plasma Spray external coating. The FIN II cup features three peripheral fins to improve primary stability and to prevent rotation and three holes in the proximal zone to allow for placement of 6.5mm screws in various lengths, with axial deviation up to 15 degrees.

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## FIN MB CUP

EN

Uncemented cup in Ti6Al4V grade 5 ELI (ISO 5832-3), with a Titanium Y 367 external coating.


The FIN MB cup features the same morphological characteristics as the FIN II cup except for the three cranial fins. Primary and secondary stability, however, is maintained by the cup's geometry and porous coating.

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## FIN DMD CUP

EN

Uncemented primary cup, made of Titanium Ti6Al4V ELI (ASTM F3001) by Additive Manufacturing. The cup is characterized by a lattice structure that replicates that of a diamond. The size and morphology of the cells provide an ideal habitat for the osteoconduction process, thus accelerating osseointegration.

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## FIN DMD MULTIHOLE CUP

EN

An uncemented revision cup made of Titanium Ti6Al4V ELI (F3001) by Additive Manufacturing, the FIN DMD MULTIHOLE CUP features holes for using cancellous screws (10 holes in sizes 46, 48 and 50 and 12 holes in sizes from 52 to 72). The cup is equipped with a wide range of wedges allowing the surgeon to deal with most acetabular revision cases, even in the presence of cavitary and segmental bone defects (for further information, please, refer to FIN DMD WEDGE surgical technique: WEDGE-OT-000).

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In the case of a DMD Multihole cup implant with a Wedge, please refer to the FIN DMD WEDGE Surgical Technique (WEDGE-OT-000). EN

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### Pre-operative assessment

EN

Pre-operative assessment of the size and positioning of the acetabular component will facilitate the definition of the direction of reaming. A bilateral X-ray of the pelvis will help gauge the length of the limb, which has to be determined pre-operatively and then corrected. The pre-operative templates of the FIN SYSTEM are provided in a 15% enlarged format.

- Using a radiographic reference may help determine the enlargement of the patient's X-ray.
- Choosing a size cup by pre-operative assessment is an approximate method. The definitive size will be determined intra-operatively.

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### Preparing the acetabulum

EN

Ream the acetabulum starting with the smallest reamer suitable for the cavity. Gradually widen the acetabulum by incrementally reaming the articular cartilage until reaching the cancellous bone.

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### Measuring with the trial cup

EN

Thread the Impactor (110380420) on to the Trial Cup (1103816XXN) as previously determined by the last acetabular reamer used. Assemble the desired adaptor for the Impactor (1103803") to the impactor (110380420). Thread the trial cup of the desired size to the impactor. The trial cup has six slots to verify that the housing of the cup has been correctly prepared and is congruous with respect to the bone. If necessary, use the trial inserts directly inside the trial cup and proceed as per (7) below.

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### Assemble the impactor with chosen cup size

EN

Insert the Adaptor for Impactor (1103803") of the desired size in the Impactor (110380420). Once assembled, attach the cup to be implanted to the impactor.

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## Inserting the final cup

EN

Place the final cup in the acetabulum and impact until stable fixation is achieved. The cup should be positioned close to 45°.

NB: In the case of a Fin II cup implant, position the 3 fins in the upper quadrant (polar) of the acetabulum.

The removable rods M4 (110380426) and M5 (110380427) on the impactor aid orientation: the 45° rod must stay perpendicular to the horizontal plane, the 12° one parallel to the acetabulum.

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## Inserting the trial insert

EN

Having removed the impactor, manually place the Trial Insert (110382XXXX) in the cup's conical cavity and secure it with the assembled screw and screwdriver(110381680). Do not impact.

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EN

Having removed the trial insert, screw the Apex Plug into the central hole in the cup using the Screwdriver (110381680).

The screw holes in the cup are closed by screw hole plugs, which can be removed when necessary.

NB: the apex plugs and screw hole plugs are supplied as follows:

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## Positioning the screws

EN

Place the Drill Guide (110388701) in the screw seat and drill with the  $\varnothing$  3.2 mm Bit (L. 56 mm, 110388772 and L. 70 mm 110388771) previously mounted on the Flexible shaft (110388760). If requested by the client, the  $\varnothing$  4.5 mm Drill Bit (L. 56 mm 110388773) is also provided.

Determine the length of the screws with Depth Gauge (110388730), pick up the chosen screw with Screw Holder Forceps (110381681) and use UJ screwdriver (170010112) to insert the screw.

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## Positioning the final insert

EN

Carefully clean the taper inside the metal cup, making sure there are no particles of soft tissue which may get between cup and insert. To place the final insert (ceramic or polyethylene), place the suction cup (110380419) on the suction cup holder (110380418). Mount on conical impactor (110380420).

Once the insert is touching the taper, tap it lightly with the Insert Impactor (110388630) using the 28 mm, 32 mm or 36 mm heads (110388635; 110388636; 100388637), depending on the chosen insert.

When positioning a polyethylene liner with an anti-luxation shoulder, match the reference on the edge of the shoulder with the marking on the body of the impactor, in order to obtain the orientation chosen with the trial.

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