# DUALIS SYSTEM

COTILE A DOPPIA MOBILITÀ  
DOUBLE MOBILITY CUP  
COTILO DE DOBLE MOVILIDAD

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## DUALIS

EN

The concept of double mobility (head-insert and insert-cup) is a tried and tested means of reducing the risk of dislocation and consequent treatment. The principle is that of combining a high stability of the joint afforded by a large diameter polyethylene insert and reduction of shear forces thanks to the "low-friction" of the head-insert coupling. The Dualis cup makes it possible to combine double mobility with press-fit fixation.

## Indications

EN

Dualis is a dual mobility cup particularly indicated in patients with high risk of dislocation as neurological patients (Alzheimer's disease, stroke Outcomes cerebral, Dementia, Psychiatric Patients, Parkinson's disease), or hardly collaborating. The Dualis cup is characterized by a cemented version, indicated for cemented fixation, and cementless version, indicated for cementless implant.

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## Contraindications

EN

The hip joint surgery is specifically contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis. It is relatively contraindicated in case of: osteoporosis; patient being uncooperative or suffering from neurological disorders, unable to follow directions; systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; osteomalacia; active infection or suspected latent infection in the hip joint; distant focus of infection that could spread to the implant site; persistent acute or chronic osteomyelitis; vascular insufficiency, muscular atrophy, neuromuscular diseases; incomplete or insufficient presence of soft tissue around the hip joint; obesity; inadequate bone stock for the prostheses support or fixation; skeletal immaturity; local or disseminated neoplastic diseases; non treatable severe deformities.

## Materials

EN

Metal component:

- Uncemented DUALIS: high nitrogen content steel to ISO 5832-9. Double Ti SPS + HA OSPROVIT plasma spray coating for greater secondary stability.
- Cemented DUALIS: high nitrogen content steel to ISO 5832-9.

Insert is made of high cross-linked polyethylene (XLPE - High Cross Linked Polyethylene) according to ISO 5834-2.

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IT

In the case of a CEMENTED DUALIS implant, please refer to pag. 10

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### Pre-operative assessment

EN

Pre-operative evaluation enables the surgeon to adopt the correct approach to the operation and anticipate situation that may arise during surgery.

On the patient's most recent X-Rays perform thorough pre-op planning by using the X-Ray templates. Assess sizes of cups and components, consider 15% magnification. Assessment is approximate and must always be checked and modified during surgery if necessary.

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### Preparing the acetabulum

EN

Ream the acetabulum starting with the smallest reamer (Ref. 110384042÷110384064\*), perpendicular to the acetabulum to find the right sphericity (Fig. 1).

Use successive reamers keeping the reamer holder in the same orientation as the cup to be implanted, 40° - 45° inclination and 10° - 15° anteversion (Fig. 2).

Proceed with the preparation to reach the desired size and make sure the reamer holder is at 40° - 45° with an anteversion of 10° - 15°.

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### Measuring the acetabulum

EN

Assemble the trial cup (Ref. 110433342 ÷110433364\*) with impactor (Ref. 110381640) (Fig.3A).  
The chosen size will be that of the definitive uncemented Dualis cup (Fig. 3B).

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### Assemble the impactor with chosen cup size

EN

Mount adaptor (Ref. 110433142 ÷110433164\*), corresponding to the size of the definitive cup, on impactor (Ref. 110433001) using expansion screw (Ref. 110433006) (Fig. 4).  
Insert the screw and partially tighten it leaving the four fins of the adaptor flexible.

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## Inserting the final cup

EN

Insert it in the cup (Fig.5A), align the window on the edge of the adaptor (Ref. 110433142÷110433164\*) with the laser mark on the edge of the cup and then screw in clockwise until tight (Fig.5B).

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EN

Partially unscrew the impactor and gently tap the top metallic part of the handle to release the adaptor from the cup (Fig.5C).

If the cup has to be impacted again, impaction device (Ref. 110433003) mounted on impactor (Ref. 110381640) may be used (Fig.5D).

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### Trial Insert

EN

Choose the identified trial head (Ref. 22.2mm: 110381020\*÷110381040\*, Ref. 28mm: 110380860÷110380890) and put it in trial insert (Ref. Diam. Int. 22.2mm: 110433242÷110433264\*, Ref. Diam int. 28mm 110433444÷110433464\*) corresponding to the implanted cup and position the assembly on the definitive femur stem (Fig. 6A).

Reduce the joint using pusher device (Ref. 110433005) (Fig. 6B).

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## Definitive implant reduction

EN

Assemble the femur head and the polyethylene insert using press (Ref. 110433007). Open the press, apply the polyethylene insert and the head and turn clockwise until the head inserts with a click (Fig. 7A).

Insert the head-insert assembly on the neck of the femur stem using pusher device (Ref. 110433005) (Fig. 7B) and reduce the joint (Fig. 7C).

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In the case of a CEMENTED DUALIS implant: steps 1,2,3 are the same of Uncemented DUALIS surgical technique. EN

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### CEMENTED DUALIS

EN

In case of a cemented Dualis cup implant, the surgical technique is as for the non cemented Dualis cup, except for the steps detailed. Instrumentation and steps remain the same for both types of cup.

STEP 1, 2, 3 ARE THE SAME AS THE SURGICAL TECHNIQUE OF THE NON CEMENTED DUALIS.

Place the impactor adaptor (Ref. 110433142÷110433164\*) corresponding to the definitive cup on the cup impactor (Ref. 110433001) WITHOUT screwing in clockwise direction.

Put and spread the necessary amount of PMMA bone cement into the acetabulum.

To ensure the appropriate cement coating, the size of the definitive Dualis cup must be one size smaller than the last acetabular broach used and its relative trial cup (Ref. 110433342 ÷ 110433364\*).

The primary fixing is obtained by using bone cement: the grooves on the outer surface of the cemented cup guarantee the perfect cement anchorage to the cup surface.

Once the definitive cup is placed in place, immediately remove the excess cement that should be interposed between the impactor adaptor and the definitive cup.

When the cement has fully catalyzed, gently remove the cup impactor with the impactor adaptor.

Be careful that the definitive cup stays in place.

PROCEED BY REFERRING TO THE SURGICAL TECHNIQUE OF THE NOT CEMENTED DUALIS CUP, FROM PHASE: 6.

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### Warnings

EN

Never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant. In particular, never use stems and femoral heads produced by different manufacturers, because coupling tapers can be non-compatible. Gruppo Bioimpianti s.r.l. will not respond if any problems arised by a similar choice. The surgeon must be aware that a damaged or scratched neck can influence stem resistance and lead to a premature neck fracture.

- Handle the product with due care to avoid accidental scratches, even during surgery. An inadequate handling of the device, pre-operative and / or intra-operative, or any damage can lead to corrosion, fatigue failure or excessive wear. Use clean gloves when handling the device.
- Use screws or accessories made of suitable material. In the case of a system composed of several elements, always use the same kit elements, in case of acetabular cups made of titanium alloy (ISO 5832/3) use, if necessary, only screws of the same material.
- Follow the surgical technique provided by Gruppo Bioimpianti.
- A complete cleaning and removal of metal particles, bone cement, as well as any other debris resulting from surgery, is necessary to minimize the wear of the articular surface.
- Do not use the device if there is any damage.
- Do not alter the system in any way.
- Inadequate choice, positioning, alignment and fixation of one or more components could result in an incorrect distribution of the stresses, with a consequent reduction of life of the prosthetic implant.
- An incorrect alignment or a plant not accurate can lead to a phenomenon of excessive wear and / or failure of the system itself.
- All implantable devices for hip replacement are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use, which affects the fatigue strength and thus increases the risk of breakage.
- Not implant a component in a patient if it had already been implanted, albeit temporarily, in another patient.
- It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability
- Patient smoking may result in a delayed healing, non-healing and / or compromised stability.
- It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components.

Implantable devices for hip replacement surgery provides the surgeons the necessary tools to reduce pain and restore the Hip joint function. While this product is generally successful in attaining these goals, it cannot be expected to withstand the activity levels and function of a normal healthy joint.

### Precautions

EN

The instruments set, realized by Gruppo Bioimpianti, was designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l. except the generic instruments present in the operating rooms (drill, oscillating blades, ) could determine negative events. Any fracture or breaking of instruments shall be reported. Surgical instruments are inevitably subject to wear with normal usage and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti. It is recommended that all instruments be regularly inspected to identify any possible damage. Do not leave in the surgical site any trial, packaging, instrument or its part.

### Possible adverse effects

EN

The possible adverse effects that may result from the use of a prosthetic system, are the classics that may result from a hip arthroplasty and include the following:

- Osteolysis (progressive bone resorption). Osteolysis can be asymptomatic and therefore routine periodic radiographic examination is essential to prevent any serious future complication
- Particles leading to increased wear rates necessitating early revision.
- Early or late deep infection which may implicate the removal of the prosthesis.
- Allergic reactions.
- Intra-operative bone fractures particularly in the presence of low bone density due to osteoporosis, bone defects resulting from previous surgeries, bone resorption, or insertion of the device.
- Damage to blood vessels or hematoma
- Dislocation, migration and / or subluxation of the prosthetic components due to trauma, incorrect attachment, loss of fixation, misalignment, improper positioning, bone resorption, soft tissue laxy, sudden movement and / or natural and / or excessive activity.
- Peri-articular calcification or ossification, with or without impediment to joint mobility.
- Undesirable variation of the length of the limb.
- Reduced Range of Movement due to an incorrect choice of the components or a inadequate placement of the implant, or femoral impingement.
- Fatigue fracture of the prosthetic components can occur as results of loss of fixation or uncomplete fixation, strenuous activity, trauma, misalignment or excessive weight.
- Corrosion at the interface between components.
- Wear and / or deformation of the articular surfaces.
- Pain.
- Post-operative bone fracture post-operative and / or post-operative pain.
- Audible sounds during motion.
- Nerve damage.
- Trochanteric avulsion or non-union as result of excess muscular tension, early weight bearing, or inadequate reattachment.
- Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialization or muscle deficiencies.

### Sterilization

EN

Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1,-2,-3 while the elements manufactured in UHMWPE are sterilized with ethylene oxide (ETO) according the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioner.

### Identification and traceability

EN

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method. To guarantee external traceability even after implantation they are provided six labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

### Information

EN

For any information or enquires about this publication or anything else, contact GRUPPO BIOCPIANTI.

### Disclaimer

EN

This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons. This document does not constitute medical advice, it does not dispense medical recommendations and it does not convey any diagnostic or therapeutic information.

Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists; however Gruppo Bioimpianti excludes any liability for improper use of informatioms.

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