KORUS

S Y S T E M

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OPERATING TECHNIQUE

TÉCNICA QUIRÚRGICA

IT | EN | ES

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Enjoy Mobility

KORUS

SYSTEM

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IT

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OPERATING TECHNIQUE

EN

In this surgical technique the guidelines for the cemented and uncemented Korus stem implant are explained through the most common approaches. However the Korus instrument set allows the surgeon to select the sequence of the procedure according to his own preferences.

Indications

EN

The indications are tied to hip pathologies that require an arthroplasty to reduce or eliminate pain and / or improve joint function.

The general guidelines are:

- Non-inflammatory degenerative joint disease such as primary or secondary osteoarthritis;
- Aseptic necrosis of the femoral head;
- Rheumatoid Arthritis;
- Post-traumatic Arthritis;
- Correction of functional deformity;
- Outcomes of fractures of the femoral neck;
- Outcomes of traumatic dislocations of the hip;
- Failures of osteotomy;
- Outcomes of arthrodesis.

The stems coated with hydroxyapatite are indicated for cementless fixation.

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Contraindications

EN

The hip joint surgery is absolutely contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis. It is relatively contraindicated in case of: Osteoporosis; Patient uncooperative or suffering from neurological disorders, unable to follow directions; Systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; Neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; Osteomalacia; Active infection or suspected latent infection in the hip joint; Distant focus of infection that could spread to the implant site; Vascular insufficiency, muscular atrophy, neuromuscular diseases; Incomplete or insufficient presence of soft tissue around the knee joint; Obesity; Inadequate bone stock for the prostheses support or fixation; Skeletal immaturity; Local or disseminated neoplastic diseases; Incorrigible severe deformities.

Materials

EN

Uncemented Korus, uncemented modular Korus, uncemented collared Korus: Titanium alloy Ti6Al4V grade 5 ELI (ISO 5832/3).
Finishing: HA coating.
Uncemented modular KORUS Necks: cobalt chrome alloy (ISO 5832/4).
Cemented Korus and cemented KORUS with distal centralizer: High nitrogen content stainless steel (ISO 5832/9).
Finishing: mirror finishing.
Distal Centralizer: PMMA and/or Polyethylene* (UHMWPE, ISO 5834/2).

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In the case of a CEMENTED KORUS STEM implant, please, refer to page n.12.

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Pre-operative planning

EN

Pre-operative evaluation enables the surgeon to adopt the correct approach to the operation and anticipate situations that may arise during surgery. The goal is to choose the size and the position of the femoral stem in the antero-posterior and lateral plane, based on elements from the patient's history, radiographic analysis and physical examinations. Choosing a stem by pre-operative evaluation is an approximate method. The final size of the stem will be determined intraoperatively. A bilateral X-ray of the pelvis will help to evaluate the leg length discrepancy, which must be corrected during surgery. To identify the leg length discrepancy radiographically, draw a reference line through the bottom of the ischium and measure the distance from the lesser trochanter landmark to the reference line. The difference is the radiographic leg length discrepancy. The leg length can be corrected by using the femoral heads, available in different lengths. Korus stem pre-operative templates are available with 15% magnification.

Use of an X-ray reference may help to determine the enlargement of the patient's X-ray. Align the template along the anatomic axis of the femur keeping the stem in the femoral channel and the centre of rotation of the femoral head correctly positioned.

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Femoral neck osteotomy

EN

Having freed the femoral neck, proceed with the osteotomy by resection approximately 1 cm above the lesser trochanter with an angle of around 45°. Perform the osteotomy, taking care to maintain the correct angle.

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Femoral canal preparation

EN

Before proceeding with the Broaches, the medial part of the great trochanter can be removed by using the Great Trochanter Chisel (Ref. L.180mm 110381011 or Ref. L.230mm 110381013).

This creates a starting cavity for Broaches, thus facilitating their insertion in the axial position.

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Broaching the femoral canal

EN

Once the Broach insertion has been carried out with the Great Trochanter Chisel, use the Solid Reamer (Ref. 110381400) to identify the femoral canal.

Then mount the Broach Handle (Ref. 120411111) on the Broach (Ref. 12041010*). Begin with the smallest Broach and progressively enlarge the metaphyseal cavity until complete stability is achieved, without reaching cortical contact in the femoral canal. Correct insertion of the first Broach is of vital importance: it must be lateral to the great trochanter and axial in the femoral canal to avoid a varus position of the implant. When the Broach can proceed no further along the femoral canal, the size used will be the size of the stem to implant.

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Optional for collared Korus: Calcar Reaming

EN

Leaving the last broach in place, start reaming with the Calcar Reamer (Ref. Diam. 45mm 120411116, Ref. Diam. 52mm 120411117). Two calcar reamer are available to help achieve the appropriate diameter.

Reaming may be done manually with the "T" handle or with a motor, using the Zimmer Adaptor.

Place the Calcar Reamer of the chosen diameter on the post of the broach.

Keep slowly reaming until a flat resection surface is achieved.

The calcar reaming should allow an optimised fit of the collar on the calcar.

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Trial Reduction

EN

Once the Broach has been determined, the last Broach used is left "in situ" and the trial reduction can be started. Insert the Trial Head Adaptor (Ref. 120411106 or 120411107) and subsequently the Trial Heads (Ref. 110380***) available in three different diameters (28mm, 32mm, 36mm), each of them available in four different lengths (short, medium, long, extra-long).

Proceed with the reduction by joining the Broach, combined with the Trial Head, to the previously implanted Cup and check the stability through the full range of motion, the leg length and the tension of the soft tissue.

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Stem Insertion

EN

Once the trial reduction has been executed, remove the Trial Head, the Trial Head Adaptor and then the Broach.

Proceed with the stem insertion using the Stem Impactor (Ref. 12041108) which allows the correct positioning.

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Positioning the femoral head

EN

Having inserted the stem, it is advisable to repeat the trial reduction (see point 5) with the Trial Head to choose the size of the definitive head. If the femoral head is metallic, use the Head Impactor (Ref. 110380800) combined with the Head Impactor Adaptor (Ref. 1103808**) chosen on the basis of the head's diameter. If the femoral head is ceramic, proceed with a manual insertion and make a half clockwise turn.

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Stem removal

EN

If Korus stem needs to be removed, use the Stem Extractor Impactor (Ref. 110380820).

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CEMENTED KORUS STEM implant: SECTIONS 1, 2 AND 3 ARE THE SAME AS THE SURGICAL TECHNIQUE FOR THE NON-CEMENTED KORUS STEM.

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CEMENTED KORUS STEM

EN

In the case of a cemented Korus stem implant the surgical technique is as for the non-cemented Korus stem except for the actions detailed below.

Once the Broach insertion has been carried out with the Great Trochanter Chisel, use the Solid Reamer (Ref. 110381400) to identify the endo-medullary canal.

Then mount the Broach Handle (Ref. 120411111) on the Broach (Ref. 12041010**). Begin with the smallest Broach and progressively enlarge the cavity until stability is achieved, without making contact with the cortical. Correct insertion of the first Broach is of vital importance: it must be lateral to the great trochanter and axial in the femoral canal to avoid a varus position of the implant.

When the Broach can proceed no further along the femoral canal, the last Broach used is left "in situ" and the trial reduction can be started.

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Sizing

EN

Once the trial reduction has been completed, remove (in the following order) the Trial Head, the Trial Head Adaptor and lastly the Broach. After putting the necessary amount of cement in the medullary channel, insert the definitive stem chosen up to the black mark on its surface.

The size of the definitive stem to insert will be one size smaller than the last Broach used.

(For example, if the last Broach used is size 6, the stem to implant will be size 5).

Proceed as per the surgical technique for the non-cemented Korus stem.

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DISTAL CENTRALIZER

EN

In the case of a cemented Korus with distal centralizer stem implant the surgical technique is as for the cemented Korus stem.

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
CEMENTED KORUS WITH DISTAL CENTRALIZER

REF	SIZE	CCD (A)	NECK LENGTH (B)	STEM LENGTH (C)	STEM LENGTH (E)	OFFSET (D)
120430021	1	135°	41mm	107.7mm	127mm	37mm
120430022	2	135°	41mm	118mm	137mm	38mm
120430023	3	135°	41mm	122.7mm	142mm	38.5mm
120430024	4	135°	41mm	128mm	147mm	39.5mm
120430025	5	135°	41mm	132.9mm	152mm	40mm
120430026	6	135°	41mm	137.7mm	157mm	40.5mm
120430027	7	135°	41mm	143mm	162mm	41.5mm
120430028	8	135°	41mm	148mm	167mm	42mm
120430029*	9	135°	41mm	157.6mm	177mm	43mm

MODULAR KORUS

REF	SIZE	STEM LENGTH (E)
120410602*	2	137mm
120410603*	3	142mm
120410604*	4	147mm
120410605*	5	152mm
120410606*	6	157mm
120410607*	7	162mm
120410608*	8	167mm
120410609*	9	177mm

MODULAR KORUS - NECKS

REF	SIZE	TYPE
120420001*	SHORT	STANDARD
120420002*	LONG	STANDARD
120420011*	SHORT	VARUS/VALGUS 8°
120420012*	LONG	VARUS/VALGUS 8°
120420021*	SHORT	VARUS/VALGUS 15°
120420022*	LONG	VARUS/VALGUS 15°
120420031*	SHORT	ANTE/RETRO 8°
120420032*	LONG	ANTE/RETRO 8°
120420041*	SHORT	ANTE/RETRO 14°
120420042*	LONG	ANTE/RETRO 14°

CEMENTED KORUS WITH DISTAL CENTRALIZER

REF	SIZE	CCD (A)	NECK LENGTH (B)	STEM LENGTH (C)	STEM LENGTH (E)	OFFSET (D)
120430031	1	125°	43.5mm	108.3mm	127mm	44mm
120430032	2	125°	42.1mm	119.6mm	137mm	45mm
120430033	3	125°	42mm	124.5mm	142mm	45.5mm
120430034	4	125°	41.8mm	130mm	147mm	46.5mm
120430035	5	125°	43.6mm	133mm	152mm	47mm
120430036	6	125°	41.7mm	140mm	157mm	47.5mm
120430037	7	125°	41.6mm	145.3mm	162mm	48.5mm
120430038	8	125°	41.5mm	150mm	167mm	49mm
120430039*	9	125°	43.9mm	158mm	177mm	50mm

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


Warnings EN

Never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant. In particular, never use stems and femoral heads produced by different manufacturers, because coupling tapers can be non-compatible. Gruppo Biomeptants S.r.l. will not respond for any problems arisen by a similar choice. The surgeon must be aware that a damaged or scratched neck can influence stem resistance and lead to a premature neck fracture.
Handle the product with due care to avoid accidental scratches, even during surgery. An inadequate handling of the device, pre-operative and / or intra-operative, or any damage can lead to corrosion, fatigue failure or excessive wear. Use clean gloves when handling the device. Use screws or accessories made of suitable material; in the case of a system composed of several elements, always use the same kit elements. In case of acetabular cups made of titanium alloy (ISO 5832/3) use, if necessary, only screws of the same material. Follow the surgical technique provided by Gruppo Biomeptants. A complete cleaning and removal of metal particles, bone cement, as well as any other debris resulting from surgery, it is necessary to minimize the wear of the articular surface. Do not use the device if there is any damage. Do not alter the system in any way. Inadequate choice, positioning, alignment and fixation of one or more components could result in an incorrect distribution of the stresses, with a consequent reduction of life of the prosthetic implant. An incorrect alignment or a plant not accurate can lead to a phenomenon of excessive wear and / or failure of the system itself. All implantable devices for hip replacement are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use which affects the fatigue strength and thus increases the risk of breakage. Do not implant a component in a patient if it had already been implanted, albeit temporarily, in another patient. Use ALWAYS the appropriate locking screw for the SMR stem (revision stem) and MRP (humoral stem). In the revision stem SMR, the locking screw must have the same length of the proximal component. In the MRP stem, the locking screw must have the same length of the proximal component. As regards the acetabular wedges of FIN DMD system fit the wedge/cup interface using the cemented fixation and fit the wedge/remaining acetabular bone using the cementless fixation by press-fit and 6.5mm screws, which must be passed through the holes present on the same wedge. Never use these screws between the two components (through cup and wedge); the wedge need to be fixed to the cup only using bone cement, without the use of the screws. Always use the two pegs and the concious screw with the Mails Iripod cup. Make sure that the two pegs of the Duais Iripod cup are correctly placed in their seat, and that they are not protruding toward the inner surface of the cup, to avoid damage of the mobile ALP insert. In case it is necessary to extract and re-position the CCo mer for the Dual Mobility of the FIN System, never reuse the extracted component, but get a new one. It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability. Patient smoking may result in a delayed healing, non-healing and/or compromised stability. It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components. The surgeons and other users need to be aware that, if explanted, the retrieved devices must be treated as contaminated and containing products and that they need to be disposed according to the specific instructions of the manufacturer. The devices retrieved may result in patient injury. Implantable devices for hip replacement surgery provides the surgeons the necessary tools to reduce pain and restore the hip joint function. While this product is generally successful in attaining these goals, it cannot be expected to withstand the activity levels and function of a normal healthy joint.

Possible adverse effects EN

The possible adverse effects that may result from the use of a prosthetic system, are the classics that may result from a hip arthroplasty and include the following: Osteolysis (progressive bone resorption). Osteolysis can be asymmetric and therefore routine periodic radiologic examination is essential to prevent any serious future complication. Particles leading to increased wear rates necessitating early revision. Early or late deep infection which may necessitate the removal of the prosthesis. Allergic reactions: Intra-operative fractures of the bone. Intra-operative fractures in particular in the presence of low bone density due to osteoporosis, bone defects resulting from previous surgeries, bone resorption, or insertion of the device; Damage to blood vessels or hematoma. Dislocation, migration and / or sublucazione of the prosthetic components due to trauma, incorrect attachment, loss of fixation, misalignment, improper positioning, bone resorption, soft tissue early sudden movement and / or natural and / or excessive activity; Peri-articular calcification or ossification, with or without impediment to joint mobility; Undesirable variation of the length of the limb; Reduced Range of Movement due to an incorrect choice of the components or a inadequate placement of the implant; or failure of movement; Fatigue fracture of the prosthetic components can occur as results of loss of fixation or incomplete fixation, stress reaction, activity, trauma, misalignment or excessive weight. Corrosion at the interface between components. Wear and / or deformation of the articular surfaces. Pain; Post-operative bone fracture post-operative and / or post-operative pain; Audible sounds during motion; Nerve damage; Trochanteric avulsion or non-union as result of excess muscular tension, early weight bearing, or inadequate attachment; Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialization or muscle deficiencies.

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Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1, -2, -3 while the elements manufactured in UHMWPE are sterilized with ethylene oxide (ETO) according to the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioning.

Surgical Instruments	EN
----------------------	----

The instruments sets related to femoral and acetabular components realized by Gruppo Bioimpianti are designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. Such instruments are supplied either as part of an operating set or dedicated trays or individually and are identified by a code and a lot number, reported on the label when present and on the device if possible. Surgical instruments must be kept in hygiene-controlled premises, if possible in dedicated conservation units (surgical instrument trays). The materials used for manufacturing are biocompatible and guarantee good resistance to corrosion if stored in appropriate conditions. It should be remembered that organic fluids can cause local corrosion, so it is good practice to clean instruments after use, removing traces of blood or other organic fluids which might prove aggressive to the material if left in such conditions. The instruments, whether part of an operating set or supplied individually, are supplied non-sterile. Gruppo Bioimpianti S.r.l. performs on each instrument a washing cycle according to the applicable operative instructions. Gruppo Bioimpianti underlines that it is responsibility of the hospital facility to perform on each instrument, before each surgery, the following sequence of processes, that must be performed appropriately according to the internal procedures approved by each facility: decontamination, washing (manual or ultrasonic), re-washing, drying, control of absence of residual dirt or other impurities, control of lumens/blades, control of the integrity of the instruments and sterilization according to a sterilization process conforming to the relevant standards (EN 556, ISO 11137-1, -2, -3 and ISO 11135-1 and UNI EN ISO 17665-1). Gruppo Bioimpianti S.r.l. manufactures instruments that can be sterilized several times, even using the simple systems generally found in hospitals (steam sterilizers). Should the instrument be impaired or show signs of superficial alteration, it should not be used. Any fracture or breaking of instruments shall be reported to the company. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades ...) could determine negative events. The instruments, if not carefully operated during the surgery, could cause damages to the implants. Surgical instruments are inevitably subject to wear with normal use and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti S.r.l. It is recommended that all instruments be regularly inspected to identify any possible damage. WARNING: if the instruments are intended to be connected with an active medical device (electrical engine), the compatibility between two elements of the connection must be verified. Do not leave in the surgical site any trial component, packaging, instrument or part of it. Being medical devices, the instruments must be strictly used by medical personnel, which must be sufficiently familiar with their use before performing the surgery.

Identification and traceability	EN
---------------------------------	----

All products manufactured by Gruppo Bioimpianti S.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method, to guarantee external traceability even after implantation they are provided six labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

Information	EN
-------------	----

For any information or enquires about this publication or anything else, contact GRUPPO BIOIMPIANTI.

Disclaimer	EN
------------	----

This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons. This document does not constitute medical advice, it does not dispense medical recommendations and it does not cover any diagnostic or therapeutic information. Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists; however Gruppo Bioimpianti excludes any liability for improper use of informations.

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